Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 1/18	

Trunk support

Trunk support narrow back

If the trunk support needs to be mounted on a seat with a narrow back, assembly set 103425-99-0 will be required. The assembly procedure is shown on page 6 of these instructions.

List of parts that should be in the set

Always included:

2 pcs Brackets, right and left.

Depending on model, the following will be included:

2 pcs Cushions, A.
 2 pcs Covers for cushion A.
 2 pcs Cushions, E.
 2 pcs Covers for cushion E.
 2 pcs Cushions, F.
 2 pcs Covers for cushion F.

Read the instructions carefully *before* starting installation

Fig. 1. Cushion A. The set with order number 103080-99-0 contains cushions and brackets, Fig. 6. Order number 103050-99-0 also contains a cushion cover.

Fig. 2. Cushion E. The set with order number 103084-99-0 contains cushions and brackets, Fig. 6. Order number 103062-99-0 also contains a cushion cover.

Fig. 3. Cushion F. The set with order number 103082-99-0 contains cushions and brackets, Fig. 6. Order number 103056-99-0 also contains a cushion cover.

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 2/18	

Assembly details

The trunk support can be mounted at several different heights. The part of the bracket mounted on the back rest has sliding grooves permitting adjustment for height. If this is insufficient, the forks on which the cushions are mounted can be reversed by undoing the outer joint and mounting them from the opposite side of the bracket arm.

High or low mounting positions, see Figs. 4 and 5, can be achieved by dismounting the joint/locking part, Fig. 6, and switching between right and left brackets. They are then turned 180° and mounted.


Fig. 6. Switching sides on the joint/locking part and then turning them by 180°, permits high or low mounting position.


Fig. 4. Trunk support brackets in high position.


Fig. 5. Trunk support brackets in low position.

Assembly of cushion and trunk support bracket

The assembly method is the same regardless of the cushion alternative selected.

The cushion can be identified by the letter of the alphabet moulded into it. Figure 8.

The assembly kit includes a bracket for left or right-hand fixing. Figures 4 and 5.

The trunk support bracket has a locking screw on each of the fork pins. See figures 4 and 5. Unscrew the two locking screws.

Some of the cushions have a protective plastic plug in the grooves, which must be removed before mounting the cushion.

Slide the fork pins into the grooves, Fig. 7. How far in will depend on the distance desired compared to the back rest. Fix the cushion by screwing in both the locking screws.


Fig. 7. Assembly of cushion and bracket. Slide the bracket forks into the grooves on the cushion.

Fig. 8. The moulded identification letter.

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 3/18	

Attachment of accessory on Corpus

Mount the brackets in the pre-drilled holes in the back rest, Fig. 9. If the holes are missing, drill them according to the dimensions in Fig. 9. If the seat has a sliding back, drill two holes only, as shown on Fig. 10. Drill the holes with a 6.5 mm bit. Fig. 11 shows trunk supports mounted on the back.


Fig. 9. If there are no pre-drilled holes, drill them according to the given dimensions. Dimensions are given in millimetres and are measured from each edge.


Fig. 10. If the chair has a sliding back, drill the holes according to the dimensions given above. Dimensions are given in millimetres and are measured from each edge.

Fig. 11. Trunk support mounted on normal-width back.

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21		Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 4/18

Attachment of accessory on PS Jr & TS/TS Jr.

1. Adapt the Trunk supports by turning the Trunk Support arm 180°. The arm is mounted with one screw. Mount according to fig. 12.


Fig. 12. Mounting of Trunk Support Arm.

2. **PS Jr:** Mount the Trunk support on the Back Rest using the accompanying screws(M6x25). See fig. 13.


Fig. 13. Mounting the Trunk support on PS Jr.

TS/TS Jr: Dismount the bracket used when fitting on PS Jr. Dismount and remount the release button on the other side of the hinge. Remove the rubber protection on the back rest frames lower edge and fit the nuts into the groove of the back rest frame. See fig.

Mount the Trunk support on the Back Rest using the accompanying screws(M6x16). See fig. 14.

Refit the rubber protection on the back rest frames lower edge.


Fig. 14. Mounting the Trunk Support on TS/TS Jr.

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21		Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 5/18

Adjustment of accessory

Apart from adjustments to the height of the trunk support, as discussed earlier, the distance to the back rest can also be adjusted, with adaptation to individual trunk measurements.

The distance to the back support is adjusted by loosening the screws holding the cushion. Set the desired distance and fix the trunk support with the screws. See fig. 7.

There are two ways of adapting the supports to individual trunk measurements. This is made possible by adjusting three joints. Undo the joint on the joint/locking part, Fig. 6, and adjust the desired distance between the cushions. Adjust the angle with the body using the outermost joint on the trunk support arm. Tighten the screws after adjustment.

The second type of adaptation is shown in Fig. 15. There are two joints next to where the cushion is fixed to the trunk support. These allow adjustment of both the gap between the cushions and the distance from the back support.

Fig. 15. Adjusting the trunk support to the user's requirements. Adjustment using the rear joint.

Fig. 16 Preliminary adjustment of the trunk support using the outer joints.

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 6/18	

Trunk support on Corpus with narrow back

Order no.: 103425-99-0

List of parts that should be in the set

- 2 pcs Brackets, right and left.
- 2 pcs Screws M5x12.
- 4 pcs Screws M6x16.
- 2 pcs Stop.

Fitting

If the trunk support needs to be mounted on a Corpus with a narrow back, a modified bracket, order no. 103425—99-0 will be required.

1. Remove the ordinary bracket from the Trunk support, it's fitted with two screws. See fig. 17.
2. Remove the release button. See fig. 17.


Fig. 17. Remove the bracket and the release button.

3. Fit the trunk support on the modified bracket, fit the release button on the inner side of the Trunk support. See fig. 18.


Fig. 18. Fitting the modified bracket and the release button.

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 7/18	

- Remove the backrest by loosening the locking peg and the locking knobs on the back of the backrest. Push the upper part of the back support slightly forward to release it from the locking peg and lift the backrest upwards to remove it. See fig. 19.


Fig. 19. Release the Locking Peg and Knobs.

- The semicircular opening on the lower part of the bracket must straddle the shaft which holds the armrests. Screw the bracket into the metal tube supporting the plastic cover using an M5x12 screw
- Refit the back rest. See Fig. 19.

Fig. 20. Mounting the trunk support.

Adjustment

Adjustment of the accessory is performed as described earlier. See page 5.

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 8/18	

Trunk support on VS with narrow back*Order no.: 1825003*

If the trunk support needs to be mounted on a VS with a narrow back, a modified bracket, order no. 1825003 will be required.

1. On the front side of the back rest, measure and mark where to drill the holes. See figures below.
2. Drill holes $\varnothing$ 6,5mm.
3. Fit the Trunk support on the back rest using the separate delivered bracket.


Fig. 21. The Trunk Support is mounted with a separate delivered bracket (Order. no. 1825003).


Fig 22. The upper hole should be positioned 44mm. from the edge of the Back Rest. $\varnothing$ 6,5mm.


Fig. 23. The lower hole should be positioned 48mm. from the edge of the Back Rest. $\varnothing$ 6,5mm.


Fig. 24. The position of the holes, measured from the lower edge of the Back Rest.

Fig. 25. The Trunk Support fitted.

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 9/18	

Trunk support on VS with extra short back

Order no.: 1825003

If the trunk support needs to be mounted on a VS with a extra short back, a modified bracket, order no. 1825003 will be required.

1. On the front side of the back rest, measure and mark where to drill the holes. See figures below.
2. Drill holes $\varnothing$ 6,5mm.
3. Fit the Trunk support on the back rest using the separate delivered bracket.

A and B measure, at different backrest widths.

Item (w×h)	A (mm)	B (mm)
311819 (346x475)	44	48
321315 (400x475)	48	52
321323 (450x475)	48	52


Fig 27. The upper $\varnothing$ 6.5 mm. hole drilled according to the table (A mm.) from back-rest edge.


Fig 28. The lower $\varnothing$ 6.5 mm. hole drilled according to the table (B mm.) from back-rest edge.


Fig. 29. The position of the holes, measured from the lower edge of the Back Rest.

Fig. 30. The Trunk Support fitted.


Fig. 26. The Trunk Support is mounted with a separate delivered bracket (Order. no. 1825003).

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 10/18	

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Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21		Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 12/18

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
Assembly Instruction


Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 13/18	

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
Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 14/18	

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Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 15/18	

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Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 16/18	

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Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 17/18	

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Bild 25. Bålstödet monterat.

Assembly Instruction

Document No. 232672	Revision C					
Production state SERIES						
Preform Date 2008-01-21	Preform By BTA	Approved Date 9/13/2013	Approved By GHT		Page 18/18	

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